

# TM Widely-Targeted Metabolomics

A combination of an untargeted assay (with QTOF and DDA) followed by an additional run with QQQ and MRM for better annotation and accurate relative quantitation of metabolites. We typically identify/ annotate over 1,400 metabolites with an accurate relative quantitation in a variety of samples.



**Large Curated Database**  
Over 3000 metabolites are from in-house purified chemical standards



**Rigorous Quality Control**  
10 QC Indicators




**Accurate Identification**  
AB SCIEX QTOF 6600+




**Precise Quantitation**  
AB SCIEX Triple Quad 6500+

## 7 Highlights




## High Stability

CV of six internal standards over 8 projects.



## Project Experience



Number of metabolites detected from various samples, serum and plasma and fecal samples.



**Contact Us**  
[support-global@metwarebio.com](mailto:support-global@metwarebio.com)

+1(781)975-1541  
8A Henshaw Street, Woburn, MA 01801  
[www.metwarebio.com](http://www.metwarebio.com)